

Product Datasheet

Probable ATP-dependent RNA helicase DDX27 antibody (HRP) (orb243109)

Description	Probable ATP-dependent RNA helicase DDX27 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Probable ATP-dependent RNA helicase DDX27 protein (30-794AA)
Target	DDX27
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96GQ7
Expiration Date	12 months from date of receipt.