

Product Datasheet

Protein argonaute-2 antibody (HRP) (orb243133)

Description	Protein argonaute-2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Protein argonaute-2 protein (517-818AA)
Target	AGO2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKV8
Expiration Date	12 months from date of receipt.