

Product Datasheet

Protein-arginine deiminase type-2 antibody (orb243150)

Description Protein-arginine deiminase type-2 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen Recombinant Human Protein-arginine deiminase type-2 protein (1-665AA)

Target PADI2

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID [Q9Y2J8](#)

Dilution Range WB:1:500-1:5000, IHC-P:1:200-1:500, IF/ICC:1:50-1:200


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human CO...



Western blot analysis of Rat brain tissue...