


Product Datasheet

Protein-arginine deiminase type-2 antibody (orb243151)

Description	Protein-arginine deiminase type-2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Protein-arginine deiminase type-2 protein (1-665AA)
Target	PADI2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2J8
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3(Lane 1), H...