


Product Datasheet

Purine nucleoside phosphorylase antibody (orb243162)

Description	Purine nucleoside phosphorylase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Purine nucleoside phosphorylase protein (1-289AA)
Target	PNP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00491
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...