

Product Datasheet

Putative tRNA (cytidine(32)/guanosine(34)-2'-O)-methyltransferase antibody (Biotin) (orb243166)

Description	Putative tRNA (cytidine(32)/guanosine(34)-2'-O)-methyltransferase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Putative tRNA (cytidine(32)/guanosine(34)-2'-O)-methyltransferase protein (1-329AA)
Target	FTSJ1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UET6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.