



Product Datasheet

RAC-beta threonine-protein kinase antibody (orb243182)


Description	RAC-beta threonine-protein kinase antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human RAC-beta serine/threonine-protein kinase protein (2-481AA)
Target	AKT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31751
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Immunofluorescence analysis of A549 cell...



Western blot analysis of MCF-7 whole cel...