

Product Datasheet

Ras-related protein Rab-11A antibody (HRP) (orb243193)

Description	Ras-related protein Rab-11A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Ras-related protein Rab-11A protein (2-213AA)
Target	RAB11A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62491
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.