

Product Datasheet

Ras-specific guanine nucleotide-releasing factor RalGPS1 antibody (HRP) (orb243205)

Description	Ras-specific guanine nucleotide-releasing factor RalGPS1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Ras-specific guanine nucleotide-releasing factor RalGPS1 protein (1-305AA)
Target	RALGPS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5JS13
Expiration Date	12 months from date of receipt.