


Product Datasheet

Repressor of RNA polymerase III transcription MAF1 homolog antibody (orb243226)

Description	Repressor of RNA polymerase III transcription MAF1 homolog antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Repressor of RNA polymerase III transcription MAF1 homolog protein (1-256AA)
Target	MAF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H063
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 (Lane 1), ...