



## Product Datasheet

### threonine-protein kinase 24 antibody (orb243256)


|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | threonine-protein kinase 24 antibody                                   |
| <b>Species/Host</b>        | Rabbit                                                                 |
| <b>Reactivity</b>          | Human                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                     |
| <b>Immunogen</b>           | Recombinant Human Serine/threonine-protein kinase 24 protein (1-431AA) |
| <b>Target</b>              | STK24                                                                  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol            |
| <b>Form/Appearance</b>     | Liquid                                                                 |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.          |
| <b>Note</b>                | For research use only                                                  |
| <b>Isotype</b>             | IgG                                                                    |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Uniprot ID</b>          | <b>Q9Y6E0</b>                                                          |
| <b>Dilution Range</b>      | WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Hela whole cel...



Immunofluorescence staining of Hela cell...