

Product Datasheet

Set1/Ash2 histone methyltransferase complex subunit ASH2 antibody (HRP) (orb243271)

Description	Set1/Ash2 histone methyltransferase complex subunit ASH2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Set1/Ash2 histone methyltransferase complex subunit ASH2 protein (1- 534AA)
Target	ASH2L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBL3
Expiration Date	12 months from date of receipt.