


Product Datasheet

Transferrin antibody (FITC) (orb243306)

Description	Transferrin antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Human Transferrin (Native Protein)
Target	TF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02787
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109 whole cel...