



Product Datasheet

ISG15 antibody (orb243353)

Description	ISG15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Human Ubiquitin-like protein ISG15 protein (2-157AA)
Target	ISG15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05161
Dilution Range	WB:1:1000-1:5000, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunoprecipitation analysis of HeLa who...