

Product Datasheet

ISG15 antibody (HRP) (orb243356)

Description	ISG15 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Ubiquitin-like protein ISG15 protein (2-157AA)
Target	ISG15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05161
Expiration Date	12 months from date of receipt.