



Product Datasheet

UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase 110 kDa subunit antibody (orb243357)


Description	UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase 110 kDa subunit antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human UDP-N-acetylglucosamine--peptide N-acetylglucosaminyltransferase 110 kDa subunit protein (606-1022AA)
Target	OGT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15294
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...



Immunohistochemical staining of human pl...



Immunofluorescence analysis of HeLa cell...