

Product Datasheet

UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase 110 kDa subunit antibody (Biotin) (orb243358)

Description	UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase 110 kDa subunit antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human UDP-N-acetylglucosamine--peptide N-acetylglucosaminyltransferase 110 kDa subunit protein (606-1022AA)
Target	OGT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15294
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.