

Product Datasheet

Uridine-cytidine kinase 1 antibody (FITC) (orb243363)

Description	Uridine-cytidine kinase 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Uridine-cytidine kinase 1 protein (1-201AA)
Target	UCK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HA47
Expiration Date	12 months from date of receipt.