

Product Datasheet

Vacuolar fusion protein MON1 homolog A antibody (FITC) (orb243367)

Description	Vacuolar fusion protein MON1 homolog A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Vacuolar fusion protein MON1 homolog A protein (205-555AA)
Target	MON1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VX9
Expiration Date	12 months from date of receipt.