

Product Datasheet

Very long-chain acyl-CoA synthetase antibody (orb243372)

Description	Very long-chain acyl-CoA synthetase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Very long-chain acyl-CoA synthetase protein (283-620AA)
Target	SLC27A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14975
Expiration Date	12 months from date of receipt.