



Product Datasheet

Vitamin K-dependent protein S antibody (orb243383)


Description	Vitamin K-dependent protein S antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Vitamin K-dependent protein S protein (42-676AA)
Target	PROS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07225
Dilution Range	IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...



Immunofluorescence
staining of HepG2
cel...



Immunohistochemistry
of paraffin-embedde...