

Product Datasheet

Loxosceles amazonica Phospholipase D LamSicTox-alphaIC1 antibody (orb243415)

Description	Loxosceles amazonica Phospholipase D LamSicTox- alphaIC1 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Loxosceles amazonica Phospholipase D LamSicTox-alphaIC1 protein (1- 273AA)
Target	Phospholipase D LamSicTox- alphaIC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C0JAZ9
Expiration Date	12 months from date of receipt.