

Product Datasheet

Macaca mulatta C-X-C motif chemokine 10 antibody (Biotin) (orb243424)

Description	Macaca mulatta C-X-C motif chemokine 10 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Monkey
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Macaca mulatta C-X-C motif chemokine 10 protein (22-98AA)
Target	CXCL10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8MIZ1
Expiration Date	12 months from date of receipt.