

Product Datasheet

Macaca mulatta Interleukin-10 antibody (Biotin) (orb243428)

Description	Macaca mulatta Interleukin-10 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Monkey
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Macaca mulatta Interleukin-10 protein (19-178AA)
Target	IL10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51496
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.