

Product Datasheet

AIF1 antibody (HRP) (orb243441)

Description	AIF1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Allograft inflammatory factor 1 protein (2-147AA)
Target	Cytotoxin 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O70200
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.