

Product Datasheet

Baculoviral IAP repeat-containing protein 5 antibody (orb243442)

Description	Baculoviral IAP repeat-containing protein 5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Baculoviral IAP repeat-containing protein 5 protein (1-140AA)
Target	Birc5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	070201
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.