

Product Datasheet

GNAI3 antibody (Biotin) (orb24345)

Description	GNAI3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(k) subunit alpha protein
Target	GNAI3
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08754
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.