


Product Datasheet

C-C motif chemokine 2 antibody (orb243454)

Description	C-C motif chemokine 2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Mouse C-C motif chemokine 2 protein (24-96AA)
Target	Ccl2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10148
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...