

Product Datasheet

Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 antibody (Biotin) (orb243463)

Description	Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 protein (426-657AA)
Target	Samhd1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q60710
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.