


Product Datasheet

Dual specificity mitogen-activated protein kinase kinase 4 antibody (orb243466)

Description	Dual specificity mitogen-activated protein kinase kinase 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Dual specificity mitogen-activated protein kinase kinase 4 protein (2-397AA)
Target	Map2k4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P47809
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG-2(Lane 1),...