


Product Datasheet

Fas apoptotic inhibitory molecule 3 antibody (orb243474)

Description	Fas apoptotic inhibitory molecule 3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Fas apoptotic inhibitory molecule 3 protein (18-262AA)
Target	Fcmr
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A1KXC4
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat spleen & ki...