

Product Datasheet

Guanylate cyclase activator 2B antibody (Biotin) (orb243488)

Description	Guanylate cyclase activator 2B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Guanylate cyclase activator 2B protein (22-106AA)
Target	Guca2b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O09051
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.