

Product Datasheet

GNG11 antibody (Biotin) (orb24349)

Description	GNG11 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11 protein (1-70AA)
Target	GNG11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61952
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.