


Product Datasheet

Heparin-binding growth factor 2 antibody (orb243495)

Description	Heparin-binding growth factor 2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Fibroblast growth factor 2 protein (10-154AA)
Target	Fgf2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15655
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...