

Product Datasheet

Isocitrate dehydrogenase [NADP] cytoplasmic antibody (FITC) [Out of stock] (orb243508)

Description	Isocitrate dehydrogenase [NADP] cytoplasmic antibody (FITC) [Out of stock]
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Isocitrate dehydrogenase [NADP] cytoplasmic protein (1-240AA)
Target	Idh1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O88844
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.