


Product Datasheet

ACTA1 antibody (orb243542)

Description	ACTA1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Mouse Actin, alpha skeletal muscle (Native Protein)
Target	Acta1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68134
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse brain(Lan...