

Product Datasheet

Mycobacterium tuberculosis Cholesterol oxidase antibody (Biotin) (orb243555)

Description	Mycobacterium tuberculosis Cholesterol oxidase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis Cholesterol oxidase protein (1-578AA)
Target	choD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WMV8
Expiration Date	12 months from date of receipt.