

Product Datasheet

Mycoplasma hyopneumoniae 46 kDa surface antibody (Biotin) (orb243573)

Description	Mycoplasma hyopneumoniae 46 kDa surface antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA)
Target	p46
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0J8
Expiration Date	12 months from date of receipt.