

Product Datasheet

Onchocerca volvulus OV-16 antibody (orb243578)

Description	Onchocerca volvulus OV-16 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Onchocerca volvulus OV-16 antigen protein (17-197AA)
Target	OV16
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31729
Expiration Date	12 months from date of receipt.