

Product Datasheet

Peanut Arachin Ahy-3 antibody (Biotin) (orb243583)

Description	Peanut Arachin Ahy-3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Peanut Arachin Ahy-3 protein (299-478AA)
Target	Arachin Ahy-3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q647H2
Expiration Date	12 months from date of receipt.