

Product Datasheet

Pigeon Albumin antibody (orb243594)

Description	Pigeon Albumin antibody
Species/Host	Rabbit
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Serum Albumin purified from Pigeon serum (Native Protein)
Target	α 1-acidglycoprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	I6QYU0
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.