

Product Datasheet

CHCHD8 antibody (orb2436)

Description	CHCHD8 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human CHCHD8 (21-87/87 aa)
Target	CHCHD8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10 kDa
Uniprot ID	Q9NYJ1
Expiration Date	12 months from date of receipt.