

Product Datasheet

Plasmodium falciparum L-lactate dehydrogenase antibody (Biotin) (orb243603)

Description	Plasmodium falciparum L-lactate dehydrogenase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Plasmodium falciparum L-lactate dehydrogenase protein (238-316AA)
Target	L-lactate dehydrogenase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q27743
Expiration Date	12 months from date of receipt.