

Product Datasheet

Porcine Albumin antibody (FITC) (orb243609)

Description	Porcine Albumin antibody (FITC)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	FITC
Immunogen	Serum Albumin purified from Pig serum (Native Protein)
Target	Beta-mammal toxin C _{ss} 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08835
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.