

Product Datasheet

Porphyromonas gingivalis Gingipain R1 antibody (orb243611)

Description	Porphyromonas gingivalis Gingipain R1 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Porphyromonas gingivalis Gingipain R1 protein (228-720AA)
Target	rgpA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28784
Expiration Date	12 months from date of receipt.