

Product Datasheet

Protein S100G antibody (FITC) (orb243633)

Description	Protein S100G antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Recombinant Rat Protein S100-G protein (2-79AA)
Target	S100g
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02634
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.