

Product Datasheet

GNG7 antibody (FITC) (orb24366)

Description	GNG7 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-7 protein (2-65AA)
Target	GNG7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60262
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.