

Product Datasheet

Saccharomyces cerevisiae Neutral trehalase antibody (Biotin) (orb243660)

Description	Saccharomyces cerevisiae Neutral trehalase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae Neutral trehalase protein (546-751AA)
Target	NTH1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32356
Expiration Date	12 months from date of receipt.