

Product Datasheet

Salmonella heidelberg D-serine dehydratase antibody (FITC) (orb243668)

Description	Salmonella heidelberg D-serine dehydratase antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Salmonella heidelberg D-serine dehydratase protein (1-440AA)
Target	dsdA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B4TA53
Expiration Date	12 months from date of receipt.