

Product Datasheet

Albumin antibody (Biotin) (orb243672)

Description	Albumin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Serum Albumin purified from Sheep serum (Native Protein)
Target	IgG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14639
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IP:1:200-2000, IP:1:200-2000
Expiration Date	12 months from date of receipt.